



## **Stay-put**

by **ROEKO**

Less cord – less time –  
perfect preparation line

# Less cord – less time – perfect

## Superior retraction with *Stay-put*

Time and again dentists are faced with a problem when placing a retraction cord: it is difficult to place in the sulcus and to keep it there.

*Stay-put* is the answer to this problem.

When you place *Stay-put* it stays put.

## *Stay-put* impregnated



*Stay-put* combines the advantages of an impregnated retraction cord with the adaptability of a fine metal filament. The pliable core is so effective that the cord is not only easy to place in the sulcus but it stays there. Overlapping is not necessary. Impregnation with aluminium chloride promotes quick haemostasis without causing any cardiovascular problems.

Saves chair time:

*Stay-put impregnated* stops bleeding quickly and effectively and does not need to be replaced. Quick, efficient, improved results.

# preparation line

## Advantages of *Stay-put* impregnated

- ensures quick haemostasis
- impregnation without epinephrine
- color coded for easy distinction of sizes
- shows good contrast to the gingiva
- pliable and adaptable
- can be pre-shaped
- easy to place in the sulcus
- does not lift out of the sulcus
- no overlapping required





## Applying *Stay-put* is so easy:

Cut a length of *Stay-put* corresponding to the circumference of the tooth. Overlapping is not necessary.

Place the center of the cord into the sulcus.

Then work out towards the ends of the cord, packing it into the sulcus.



*Stay-put impregnated* is available in sizes x-fine, fine, medium and thick.

## *Stay-put non impregnated*

The proven non impregnated version. Provides you with the flexibility to impregnate the cords at the time of use to suit the conditions and techniques.



# Accessories

Roeko's **Retraction Cord Packers** facilitate placing the retraction cord. Their balanced notches hold the cord in place.



**Have you already tried Comprecap?**

Use *Comprecap* in combination with *Stay-put*. *Comprecap* exerts more pressure onto the retraction cord and the gingival area and stops bleeding naturally by compression. After removal the result is a clean, dry, open sulcus and a well defined margin. Best conditions for a perfect impression.



**Comprecap anatomic**

Anatomically formed compression caps have semicircled shaped spaces on two opposite sides and can therefore be easily placed on adjacent teeth.



When using *Stay-put* and *Comprecap* you'll get best conditions for impression taking with AFFINIS. Impressively free of voids or margin distortions – perfect impressions.



# Order information

## **Stay-put impregnated**

|        |             |
|--------|-------------|
| x-fine | REF 521 020 |
| fine   | REF 521 021 |
| medium | REF 521 022 |
| thick  | REF 521 023 |

## **Stay-put non impregnated**

|        |             |
|--------|-------------|
| x-fine | REF 521 000 |
| fine   | REF 521 001 |
| medium | REF 521 002 |
| thick  | REF 521 003 |

## **Comprecap anatomic**

|              | <b>Content</b> | <b>REF</b> |
|--------------|----------------|------------|
| Set (1/3/5)  | 120 pcs.       | 531 299    |
| Single sizes |                |            |
| 1            | 120 pcs.       | 531 201    |
| 3            | 120 pcs.       | 531 203    |
| 5            | 60 pcs.        | 531 205    |

## **Comprecap**

|              | <b>Content</b> | <b>REF</b> |
|--------------|----------------|------------|
| Set (1/3/5)  | 150 pcs.       | 530 000    |
| Single sizes |                |            |
| 1            | 120 pcs.       | 530 008    |
| 1            | 600 pcs.       | 531 008    |
| 3            | 120 pcs.       | 530 010    |
| 3            | 600 pcs.       | 531 010    |
| 5            | 60 pcs.        | 530 014    |
| 5            | 400 pcs.       | 531 014    |

## **Retraction Cord Packer**

|          |             |
|----------|-------------|
| smooth   | REF 520 000 |
| serrated | REF 520 001 |

Please contact us for more information or visit our website: [www.coltenewhaledent.com](http://www.coltenewhaledent.com)

### **Coltène/Whaledent AG**

Feldwiesenstraße 20  
9450 Altstätten/Switzerland  
Tel. +41 (0)71 757 53 00  
Fax +41 (0)71 757 53 01  
[office@coltenewhaledent.ch](mailto:office@coltenewhaledent.ch)

### **Coltène/Whaledent Inc.**

235 Ascot Parkway  
Cuyahoga Falls, Ohio 44223/USA  
Tel. +1 330 916 8800  
Fax +1 330 916 7077  
[productinfo@coltenewhaledent.com](mailto:productinfo@coltenewhaledent.com)

### **Coltène/Whaledent GmbH + Co. KG**

Raiffeisenstraße 30  
89129 Langenau/Germany  
Tel. +49 (0)7345 805 0  
Fax +49 (0)7345 805 201  
[productinfo@coltenewhaledent.de](mailto:productinfo@coltenewhaledent.de)

